

AESTHETIC Skin Air cooling System

*Enjoy More
Comfortable*

IC00L

- Laser Treatment Powerful cooling air temperature -30°C
- Strong wind speed 1-10 level adjustable
- Cool skin to numbness Avoid thermal injury
- Minimize pain and ensure a better result during laser





AESTHETIC



IMPORT

Compressors imported from Japan



WIND SPEED

Super wind speed blower



NOISE REDUCTION

Independent noise reduction system




-30°C

-30°C super-low temperature refrigeration



START STOP

One-key start-stop button



TRANSFORM

Stepless speed regulation, temperature adjustment knob



After Pico Laser, Co2 or Diode Laser, use Air cooling machine reducing redness, make the treatment process more effective and smooth.

Thermal damage caused by OPT fast hair removal treatment, you can use Air cooling machine lower skin temperature.



TREATMENT HANDLE



Round Adapter

For reducing the skin temperature of small treatment area, such as eyebrow, underarm, forehead.



Middle Square Adapter

For reducing the skin temperature of middle square area, especially for hair removal treatment, such as Arm, underarm, leg.



Large Square Adapter

For reducing the skin temperature of large treatment area, such as thigh, belly, especially for hair removal treatment.



SYSTEM INTERFACE



Defrost



Water flow alarm



Power on



Air volume



Temperature



SPECIFICATIONS

**System type**

ICOOL skin cooling system

Total power

1000W

Minimum temperature

-30°C

Refrigerant

R404A

Refrigeration intensity

10 levels

Cold air duct length

250mm

Cold air duct weight

200g

Air outlet size

ø6mm 10mm 10mm*1.5mm

LCD screen size

1.54 inches

Machine size

810*510*961mm 55kg

Packaging size

860*560*1011mm 97.38kg

Working voltage

AC220V±10%50Hz

AC110V±10%60Hz

BEST PARTNER



Reduce skin pain and the risk of burns during hair removal



When removing freckles, moles, and tattoos, it can effectively reduce skin pain and reduce the risk of blistering



When removing scars, moles, and skin rejuvenation, it can accelerate the heat dissipation of the skin and reduce pain



When removing freckles, moles, and tattoos, it can effectively reduce skin pain and reduce the risk of blistering